



Due to its characteristics of pore size, surface area, percentage of covering (%C), and the kind of silica it is build of, it is the suitable alternative to Hypersil ODS packings. Its chromatographic behavior exactly reproduces the one of this popular packing, being able to transfer the chromatographic methods without any kind of adjustment.

## 5 Microns Packing

### Analytical columns 45 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	5	3	0.46	TR-011000	
Hyperpack	ODS	5	4	0.46	TR-011001	
Hyperpack	ODS	5	5	0.46	TR-011002	
Hyperpack	ODS	5	10	0.46	TR-011003	
Hyperpack	ODS	5	15	0.46	TR-011004	
Hyperpack	ODS	5	20	0.46	TR-011005	
Hyperpack	ODS	5	25	0.46	TR-011006	
Hyperpack	C8	5	3	0.46	TR-011021	
Hyperpack	C8	5	4	0.46	TR-011022	
Hyperpack	C8	5	5	0.46	TR-011023	
Hyperpack	C8	5	10	0.46	TR-011024	
Hyperpack	C8	5	15	0.46	TR-011025	
Hyperpack	C8	5	20	0.46	TR-011026	
Hyperpack	C8	5	25	0.46	TR-011027	

## 5 Microns Packing

### Analytical columns 40 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	5	3	0.4	TR-411000	
Hyperpack	ODS	5	4	0.4	TR-411001	
Hyperpack	ODS	5	5	0.4	TR-411002	
Hyperpack	ODS	5	10	0.4	TR-411003	
Hyperpack	ODS	5	15	0.4	TR-411004	
Hyperpack	ODS	5	20	0.4	TR-411005	
Hyperpack	ODS	5	25	0.4	TR-411006	
Hyperpack	C8	5	3	0.4	TR-410081	
Hyperpack	C8	5	4	0.4	TR-410082	
Hyperpack	C8	5	5	0.4	TR-410083	
Hyperpack	C8	5	10	0.4	TR-410084	
Hyperpack	C8	5	15	0.4	TR-410085	
Hyperpack	C8	5	20	0.4	TR-410086	
Hyperpack	C8	5	25	0.4	TR-410087	

## 5 Microns Packing

### Microbore columns 2.1 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	5	3	0.21	TR-010272	
Hyperpack	ODS	5	5	0.21	TR-010273	
Hyperpack	ODS	5	10	0.21	TR-010274	
Hyperpack	ODS	5	15	0.21	TR-010275	
Hyperpack	ODS	5	20	0.21	TR-010276	
Hyperpack	ODS	5	25	0.21	TR-010277	
Hyperpack	C8	5	3	0.21	TR-011028	
Hyperpack	C8	5	5	0.21	TR-011029	
Hyperpack	C8	5	10	0.21	TR-011030	
Hyperpack	C8	5	15	0.21	TR-011031	
Hyperpack	C8	5	20	0.21	TR-011032	
Hyperpack	C8	5	25	0.21	TR-011033	

## 5 Microns Packing

### Microbore columns 3.0 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	5	3	0.3	TR-010278	
Hyperpack	ODS	5	5	0.3	TR-010279	
Hyperpack	ODS	5	10	0.3	TR-010280	
Hyperpack	ODS	5	15	0.3	TR-010281	
Hyperpack	ODS	5	20	0.3	TR-010282	
Hyperpack	ODS	5	25	0.3	TR-010283	
Hyperpack	C8	5	3	0.3	TR-011160	
Hyperpack	C8	5	5	0.3	TR-011034	
Hyperpack	C8	5	10	0.3	TR-011035	
Hyperpack	C8	5	15	0.3	TR-011036	
Hyperpack	C8	5	20	0.3	TR-011037	
Hyperpack	C8	5	25	0.3	TR-011038	

## 5 Microns Packing

### Semi-preparative columns HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	5	10	0.78	TR-010284	
Hyperpack	ODS	5	15	0.78	TR-010285	
Hyperpack	ODS	5	25	0.78	TR-010286	
Hyperpack	ODS	5	10	1.00	TR-010287	
Hyperpack	ODS	5	15	1.00	TR-010288	

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	5	25	1.00	TR-010289	
Hyperpack	ODS	5	5	2.12	TR-010290	
Hyperpack	ODS	5	10	2.12	TR-010291	
Hyperpack	ODS	5	15	2.12	TR-010292	
Hyperpack	ODS	5	25	2.12	TR-010293	
Hyperpack	C8	5	10	0.78	TR-011039	
Hyperpack	C8	5	15	0.78	TR-011040	
Hyperpack	C8	5	25	0.78	TR-011041	
Hyperpack	C8	5	10	1.00	TR-011042	
Hyperpack	C8	5	15	1.00	TR-011043	
Hyperpack	C8	5	25	1.00	TR-011044	
Hyperpack	C8	5	5	2.12	TR-011045	
Hyperpack	C8	5	10	2.12	TR-011046	
Hyperpack	C8	5	15	2.12	TR-011047	
Hyperpack	C8	5	25	2.12	TR-011048	

## 5 Microns Packing Novafix™ Cartridges 4 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	5	7,5	0.4	TR-010294	
Hyperpack	ODS	5	10	0.4	TR-010295	
Hyperpack	ODS	5	15	0.4	TR-010296	
Hyperpack	ODS	5	25	0.4	TR-010297	
Hyperpack	C8	5	7,5	0.4	TR-011049	
Hyperpack	C8	5	10	0.4	TR-011050	
Hyperpack	C8	5	15	0.4	TR-011051	
Hyperpack	C8	5	25	0.4	TR-011052	

## 3 Microns Packing Ultrarapid columns 4.6 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	3	3	0.46	TR-010298	
Hyperpack	ODS	3	4	0.46	TR-010299	
Hyperpack	ODS	3	5	0.46	TR-010300	
Hyperpack	ODS	3	10	0.46	TR-010301	
Hyperpack	ODS	3	15	0.46	TR-010302	
Hyperpack	ODS	3	20	0.46	TR-010303	
Hyperpack	ODS	3	25	0.46	TR-010304	
Hyperpack	C8	3	3	0.46	TR-011053	
Hyperpack	C8	3	4	0.46	TR-011054	
Hyperpack	C8	3	5	0.46	TR-011055	
Hyperpack	C8	3	10	0.46	TR-011056	
Hyperpack	C8	3	15	0.46	TR-011057	
Hyperpack	C8	3	20	0.46	TR-011058	
Hyperpack	C8	3	25	0.46	TR-011059	

## 3 Microns Packing Analytical columns 4.0 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	3	3	0.4	TR-410298	
Hyperpack	ODS	3	4	0.4	TR-410299	
Hyperpack	ODS	3	5	0.4	TR-410300	
Hyperpack	ODS	3	10	0.4	TR-410301	
Hyperpack	ODS	3	15	0.4	TR-410302	
Hyperpack	ODS	3	20	0.4	TR-410303	
Hyperpack	ODS	3	25	0.4	TR-410304	
Hyperpack	C8	3	4	0.4	TR-011060	
Hyperpack	C8	3	5	0.4	TR-011061	
Hyperpack	C8	3	10	0.4	TR-011062	
Hyperpack	C8	3	15	0.4	TR-011063	
Hyperpack	C8	3	20	0.4	TR-011064	
Hyperpack	C8	3	25	0.4	TR-011065	

## 3 Microns Packing Microbore columns 2.1 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	3	3	0.21	TR-010305	
Hyperpack	ODS	3	5	0.21	TR-010306	
Hyperpack	ODS	3	10	0.21	TR-010307	
Hyperpack	ODS	3	15	0.21	TR-010308	
Hyperpack	ODS	3	20	0.21	TR-010309	
Hyperpack	ODS	3	25	0.21	TR-010310	
Hyperpack	C8	3	3	0.21	TR-011066	
Hyperpack	C8	3	5	0.21	TR-011067	
Hyperpack	C8	3	10	0.21	TR-011068	
Hyperpack	C8	3	15	0.21	TR-011069	
Hyperpack	C8	3	20	0.21	TR-011070	
Hyperpack	C8	3	25	0.21	TR-011071	

## 3 Microns Packing Microbore columns 3.0 mm i.d. HYPERPACK

Packing	Funct.	µm	Length		Diameter	Cat.Nbr.
			µm	cm		
Hyperpack	ODS	3	3	0.3	TR-010311	
Hyperpack	ODS	3	5	0.3	TR-010312	
Hyperpack	ODS	3	10	0.3	TR-010313	
Hyperpack	ODS	3	15	0.3	TR-010314	
Hyperpack	ODS	3	20	0.3	TR-010315	
Hyperpack	ODS	3	25	0.3	TR-010316	
Hyperpack	ODS	3	3	0.3	TR-011072	
Hyperpack	ODS	3	5	0.3	TR-011073	
Hyperpack	ODS	3	10	0.3	TR-011074	
Hyperpack	ODS	3	15	0.3	TR-011075	
Hyperpack	ODS	3	20	0.3	TR-011076	
Hyperpack	ODS	3	25	0.3	TR-011077	



It reproduces with total fidelity the chromatographic behavior of the columns Hypersil BDS C18. Available in 3 and 5µm and in all length and diameter configurations.

## 5 Microns Packing

### Analytical columns 4.6 mm i.d. HYPERPACK BASIC

Packing	Funct.	µm	Length Diameter		Cat.Nbr.
			cm	cm	
Hyperpack BASIC	ODS	5	3	0.46	TR-011007
Hyperpack BASIC	ODS	5	4	0.46	TR-011008
Hyperpack BASIC	ODS	5	5	0.46	TR-011009
Hyperpack BASIC	ODS	5	10	0.46	TR-011010
Hyperpack BASIC	ODS	5	15	0.46	TR-011011
Hyperpack BASIC	ODS	5	20	0.46	TR-011012
Hyperpack BASIC	ODS	5	25	0.46	TR-011013
Hyperpack BASIC	C8	5	3	0.46	TR-011108
Hyperpack BASIC	C8	5	4	0.46	TR-011109
Hyperpack BASIC	C8	5	5	0.46	TR-011110
Hyperpack BASIC	C8	5	10	0.46	TR-011111
Hyperpack BASIC	C8	5	15	0.46	TR-011112
Hyperpack BASIC	C8	5	20	0.46	TR-011113
Hyperpack BASIC	C8	5	25	0.46	TR-011114

## 5 Microns Packing

### Analytical columns 4.0 mm i.d. HYPERPACK BASIC

Packing	Funct.	µm	Length Diameter		Cat.Nbr.
			cm	cm	
Hyperpack BASIC	ODS	5	3	0.4	TR-411007
Hyperpack BASIC	ODS	5	4	0.4	TR-411008
Hyperpack BASIC	ODS	5	5	0.4	TR-411009
Hyperpack BASIC	ODS	5	10	0.4	TR-411010
Hyperpack BASIC	ODS	5	15	0.4	TR-411011
Hyperpack BASIC	ODS	5	20	0.4	TR-411012
Hyperpack BASIC	ODS	5	25	0.4	TR-411013

Packing	Funct.	µm	Length Diameter		Cat.Nbr.
			cm	cm	
Hyperpack BASIC	C8	5	3	0.4	TR-411108
Hyperpack BASIC	C8	5	4	0.4	TR-411109
Hyperpack BASIC	C8	5	5	0.4	TR-411110
Hyperpack BASIC	C8	5	10	0.4	TR-411111
Hyperpack BASIC	C8	5	15	0.4	TR-411112
Hyperpack BASIC	C8	5	20	0.4	TR-411113
Hyperpack BASIC	C8	5	25	0.4	TR-411114

## 5 Microns Packing

### Microbore columns 21 mm i.d. HYPERPACK BASIC

Packing	Funct.	µm	Length Diameter		Cat.Nbr.
			cm	cm	
Hyperpack BASIC	ODS	5	3	0.21	TR-010317
Hyperpack BASIC	ODS	5	5	0.21	TR-010318
Hyperpack BASIC	ODS	5	10	0.21	TR-010319
Hyperpack BASIC	ODS	5	15	0.21	TR-010320
Hyperpack BASIC	ODS	5	20	0.21	TR-010321
Hyperpack BASIC	ODS	5	25	0.21	TR-010322
Hyperpack BASIC	C8	5	3	0.21	TR-011115
Hyperpack BASIC	C8	5	5	0.21	TR-011116
Hyperpack BASIC	C8	5	10	0.21	TR-011117
Hyperpack BASIC	C8	5	15	0.21	TR-011118
Hyperpack BASIC	C8	5	20	0.21	TR-011119
Hyperpack BASIC	C8	5	25	0.21	TR-011120

## 5 Microns Packing

### Microbore columns 3.0 mm i.d. HYPERPACK BASIC

Packing	Funct.	µm	Length Diameter		Cat.Nbr.
			cm	cm	
Hyperpack BASIC	ODS	5	3	0.3	TR-010323
Hyperpack BASIC	ODS	5	5	0.3	TR-010324
Hyperpack BASIC	ODS	5	10	0.3	TR-010325
Hyperpack BASIC	ODS	5	15	0.3	TR-010326
Hyperpack BASIC	ODS	5	20	0.3	TR-010327
Hyperpack BASIC	ODS	5	25	0.3	TR-010328
Hyperpack BASIC	C8	5	3	0.3	TR-011121
Hyperpack BASIC	C8	5	5	0.3	TR-011122
Hyperpack BASIC	C8	5	10	0.3	TR-011123
Hyperpack BASIC	C8	5	15	0.3	TR-011124
Hyperpack BASIC	C8	5	20	0.3	TR-011125
Hyperpack BASIC	C8	5	25	0.3	TR-011126

## 5 Microns Packing

### Semi Preparative columns HYPERPACK BASIC

Packing	Funct.	µm	Length Diameter		Cat.Nbr.
			cm	cm	
Hyperpack BASIC	ODS	5	10	0.78	TR-010329
Hyperpack BASIC	ODS	5	15	0.78	TR-010330
Hyperpack BASIC	ODS	5	25	0.78	TR-010331
Hyperpack BASIC	ODS	5	10	1.00	TR-010332

Packing	Funct.	Length Diameter			Cat.Nbr.
		µm	cm	cm	
Hyperpack BASIC	ODS	5	15	1.00	TR-010333
Hyperpack BASIC	ODS	5	25	1.00	TR-010334
Hyperpack BASIC	ODS	5	5	2.12	TR-010335
Hyperpack BASIC	ODS	5	10	2.12	TR-010336
Hyperpack BASIC	ODS	5	15	2.12	TR-010337
Hyperpack BASIC	ODS	5	25	2.12	TR-010338
Hyperpack BASIC	C8	5	10	0.78	TR-011127
Hyperpack BASIC	C8	5	15	0.78	TR-011128
Hyperpack BASIC	C8	5	25	0.78	TR-011129
Hyperpack BASIC	C8	5	10	1.00	TR-011130
Hyperpack BASIC	C8	5	15	1.00	TR-011131
Hyperpack BASIC	C8	5	25	1.00	TR-011132
Hyperpack BASIC	C8	5	5	2.12	TR-011133
Hyperpack BASIC	C8	5	10	2.12	TR-011134
Hyperpack BASIC	C8	5	15	2.12	TR-011135
Hyperpack BASIC	C8	5	25	2.12	TR-011136

## 5 Microns Packing Cartridge System 4 mm i.d. HYPERPACK BASIC

Packing	Funct.	Length Diameter			Cat.Nbr.
		µm	cm	cm	
Hyperpack BASIC	ODS	5	7,5	0.4	TR-010339
Hyperpack BASIC	ODS	5	10	0.4	TR-010340
Hyperpack BASIC	ODS	5	15	0.4	TR-010341
Hyperpack BASIC	ODS	5	25	0.4	TR-010342
Hyperpack BASIC	C8	5	7,5	0.4	TR-011137
Hyperpack BASIC	C8	5	10	0.4	TR-011138
Hyperpack BASIC	C8	5	15	0.4	TR-011139
Hyperpack BASIC	C8	5	25	0.4	TR-011140

## 3 Microns Packing Analytical Columns 4.6 mm i.d. HYPERPACK BASIC

Packing	Funct.	Length Diameter			Cat.Nbr.
		µm	cm	cm	
Hyperpack BASIC	ODS	3	3	0.46	TR-011014
Hyperpack BASIC	ODS	3	4	0.46	TR-011015
Hyperpack BASIC	ODS	3	5	0.46	TR-011016
Hyperpack BASIC	ODS	3	10	0.46	TR-011017
Hyperpack BASIC	ODS	3	15	0.46	TR-011018
Hyperpack BASIC	ODS	3	20	0.46	TR-011019
Hyperpack BASIC	ODS	3	25	0.46	TR-011020
Hyperpack BASIC	C8	3	3	0.46	TR-011141
Hyperpack BASIC	C8	3	4	0.46	TR-011142
Hyperpack BASIC	C8	3	5	0.46	TR-011143
Hyperpack BASIC	C8	3	10	0.46	TR-011144
Hyperpack BASIC	C8	3	15	0.46	TR-011145
Hyperpack BASIC	C8	3	20	0.46	TR-011146
Hyperpack BASIC	C8	3	25	0.46	TR-011147

## 3 Microns Packing Ultrarapid Columns 4.0 mm i.d. HYPERPACK BASIC

Packing	Funct.	Length Diameter			Cat.Nbr.
		µm	cm	cm	
Hyperpack BASIC	ODS	3	3	0.4	TR-411014
Hyperpack BASIC	ODS	3	4	0.4	TR-411015
Hyperpack BASIC	ODS	3	5	0.4	TR-411016
Hyperpack BASIC	ODS	3	10	0.4	TR-411017
Hyperpack BASIC	ODS	3	15	0.4	TR-411018
Hyperpack BASIC	ODS	3	20	0.4	TR-411019
Hyperpack BASIC	ODS	3	25	0.4	TR-411020
Hyperpack BASIC	C8	3	3	0.4	TR-411141
Hyperpack BASIC	C8	3	4	0.4	TR-411142
Hyperpack BASIC	C8	3	5	0.4	TR-411143
Hyperpack BASIC	C8	3	10	0.4	TR-411144
Hyperpack BASIC	C8	3	15	0.4	TR-411145
Hyperpack BASIC	C8	3	20	0.4	TR-411146
Hyperpack BASIC	C8	3	25	0.4	TR-411147

## 3 Microns Packing Microbore Columns 2.1 mm i.d. HYPERPACK BASIC

Packing	Funct.	Length Diameter			Cat.Nbr.
		µm	cm	cm	
Hyperpack BASIC	ODS	3	3	0.21	TR-010343
Hyperpack BASIC	ODS	3	5	0.21	TR-010344
Hyperpack BASIC	ODS	3	10	0.21	TR-010345
Hyperpack BASIC	ODS	3	15	0.21	TR-010346
Hyperpack BASIC	ODS	3	20	0.21	TR-010347
Hyperpack BASIC	ODS	3	25	0.21	TR-010348
Hyperpack BASIC	C8	3	3	0.21	TR-011148
Hyperpack BASIC	C8	3	5	0.21	TR-011149
Hyperpack BASIC	C8	3	10	0.21	TR-011150
Hyperpack BASIC	C8	3	15	0.21	TR-011151
Hyperpack BASIC	C8	3	20	0.21	TR-011152
Hyperpack BASIC	C8	3	25	0.21	TR-011153

## 3 Microns Packing Microbore Columns 3.0 mm i.d. HYPERPACK BASIC

Packing	Funct.	Length Diameter			Cat.Nbr.
		µm	cm	cm	
Hyperpack BASIC	ODS	3	3	0.3	TR-010349
Hyperpack BASIC	ODS	3	5	0.3	TR-010350
Hyperpack BASIC	ODS	3	10	0.3	TR-010351
Hyperpack BASIC	ODS	3	15	0.3	TR-010352
Hyperpack BASIC	ODS	3	20	0.3	TR-010353
Hyperpack BASIC	ODS	3	25	0.3	TR-010354
Hyperpack BASIC	C8	3	3	0.3	TR-011154
Hyperpack BASIC	C8	3	5	0.3	TR-011155
Hyperpack BASIC	C8	3	10	0.3	TR-011156
Hyperpack BASIC	C8	3	15	0.3	TR-011157
Hyperpack BASIC	C8	3	20	0.3	TR-011158
Hyperpack BASIC	C8	3	25	0.3	TR-011159