





Liners for Agilent Capillary GCs

	Cat.No	Description	Similar Item No.	OD (mm)	Length (mm)	Pk
	TR-L04010	2 mm Splitless		6.5	79	1
	TR-L04010-5	2 mm Splitless		6.5	79	5
	TR-L04010-10	2 mm Splitless		6.5	79	10
	TR-L04012	2 mm Splitless, Quartz	5181-8818	6.5	79	1
	TR-L04012-5	2 mm Splitless, Quartz		6.5	79	5
	TR-L04012-10	2 mm Splitless, Quartz		6.5	79	10
	TR-L04014	4 mm Split/Splitless		6.3	79	1
	TR-L04014-5	4 mm Split/Splitless	210-3003	6.3	79	5
	TR-L04014-10	4 mm Split/Splitless		6.3	79	10
	TR-L04016	4 mm Split/Splitless, Quartz		6.3	79	1
	TR-L04016-5	4 mm Split/Splitless, Quartz		6.3	79	5
	TR-L04016-10	4 mm Split/Splitless, Quartz		6.3	79	10
	TR-L04018	4 mm Split/Splitless, with Deactivated Glass Wool	19251-60540	6.3	79	1
	TR-L04018-5	4 mm Split/Splitless, with Deactivated Glass Wool		6.3	79	5
	TR-L04018-10	4 mm Split/Splitless, with Deactivated Glass Wool		6.3	79	10
	TR-L04020	Split Liner 4 mm id with Cup	18740-80190	6.3	79	1
	TR-L04020-5	Split Liner 4 mm id with Cup		6.3	79	5
	TR-L04020-10	Split Liner 4 mm id with Cup		6.3	79	10
	TR-L04022	Split Cup Liner, 4 mm id with Deactivated Glass Wool		6.3	79	1
	TR-L04022-5	Split Cup Liner, 4 mm id with Deactivated Glass Wool		6.3	79	5
	TR-L04022-10	Split Cup Liner, 4 mm id with Deactivated Glass Wool		6.3	79	10
	TR-L04024	Split Cup Liner, 4 mm with OV-1/Chromosorb W-P 80/100	18740-60840	6.3	79	1
	TR-L04024-5	Split Cup Liner, 4 mm with OV-1/Chromosorb W-P 80/100		6.3	79	5
	TR-L04024-10	Split Cup Liner, 4 mm id with OV-1/Chromosorb W-P 80/100		6.3	79	10
	TR-L04026	Single Taper Liner, 2 mm id		6.5	79	1
	TR-L04026-5	Single Taper Liner, 2 mm id		6.5	79	5
	TR-L04026-10	Single Taper Liner, 2 mm id		6.5	79	10
	TR-L04028	Single Taper Liner, 4 mm id	5181-3316	6.5	79	1
	TR-L04028-5	Single Taper Liner, 4 mm id		6.5	79	5
	TR-L04028-10	Single Taper Liner, 4 mm id		6.5	79	10
	TR-L04030	Single Taper Liner, 4 mm id with Deactivated Glass Wool	5062-3587	6.5	79	1
	TR-L04030-5	Single Taper Liner, 4 mm id with Deactivated Glass Wool		6.5	79	5
	TR-L04030-10	Single Taper Liner, 4 mm id with Deactivated Glass Wool		6.5	79	10
	TR-L04032	Double Taper Liner, 4 mm id	5181-3315	6.5	79	1
	TR-L04032-5	Double Taper Liner, 4 mm id		6.5	79	5
	TR-L04032-10	Double Taper Liner, 4 mm id		6.5	79	10

Liners for Agilent Packed columns GCs

	Cat.No	Description	Similar Item No.	OD (mm)	Lenght (mm)	Pk
	TR-L04050	Universal Packed Inj. Liner, Deactivated	5181-3382	93		1
	TR-L04050-25	Universal Packed Inj. Liner, Deactivated		93		25
	TR-L04052	Universal Packed Inj. Liner, Not Deactivated	5080-8732	93		1
	TR-L04052-25	Universal Packed Inj. Liner, Not Deactivated		93		25







Liners for Perkin Elmer Capillary GCs

	TR-L04070	Split Liner for Autosystem	N610-1052	4	92	1
	TR-L04070-5	Split Liner for Autosystem		4	92	5
	TR-L04070-10	Split Liner for Autosystem		4	92	10
	TR-L04072	Split Liner for Autosystem with Deactivated Glass Wool		4	92	1
	TR-L04072-5	Split Liner for Autosystem with Deactivated Glass Wool		4	92	5
	TR-L04072-10	Split Liner for Autosystem with Deactivated Glass Wool		4	92	10
	TR-L04074	Split Liner for Autosystem	N-612-1372	2	92	1
	TR-L04074-5	Split Liner for Autosystem		2	92	5
	TR-L04074-10	Split Liner for Autosystem		2	92	10
	TR-L04076	Dimple Splitter	0330-5181		100	1
	TR-L04076-5	Dimple Splitter			100	5
	TR-L04076-10	Dimple Splitter			100	10
	TR-L04078	Splitless	0330-5180		100	1
	TR-L04078-5	Splitless			100	5
	TR-L04078-10	Splitless			100	10

Liners for Perkin Elmer Packed columns GCs

	Cat.No	Description	Similar Item No.	OD (mm)	Lenght (mm)	Pk
	TR-L04090	Glass Liner for Autosystem	N-610-1048		112	1
	TR-L04090-5	Glass Liner for Autosystem			112	5
	TR-L04090-10	Glass Liner for Autosystem			112	10
	TR-L04092	Glass Liner for 8000 Series, Sigma 2000	0330-2221		101	1
	TR-L04092-5	Glass Liner for 8000 Series, Sigma 2000			101	5
	TR-L04092-10	Glass Liner for 8000 Series, Sigma 2000			101	10




Liners for Thermo Finnigan Capillary GCs

	Cat.No	Description	Similar Item No.	OD (mm)	Length (mm)	Pk
	TR-L05000	Split Liner, 3 mm ID	453 20031	105	1	
	TR-L05000-5	Split Liner, 3 mm ID		105	5	
	TR-L05000-10	Split Liner, 3 mm ID		105	10	
	TR-L05002	Split Liner, 5 mm ID	453 20030	105	1	
	TR-L05002-5	Split Liner, 5 mm ID		105	5	
	TR-L05002-10	Split Liner, 5 mm ID		105	10	
	TR-L05004	Splitless Liner, 3 mm ID	453 20032	105	1	
	TR-L05004-5	Splitless Liner, 3 mm ID		105	5	
	TR-L05004-10	Splitless Liner, 3 mm ID		105	10	
	TR-L05006	Splitless Liner, 5 mm ID	453 20033	105	1	
	TR-L05006-5	Splitless Liner, 5 mm ID		105	5	
	TR-L05006-10	Splitless Liner, 5 mm ID		105	10	
	TR-L05008	For Wide Bore Columns	453 00310	105	1	
	TR-L05008-5	For Wide Bore Columns		105	5	
	TR-L05008-10	For Wide Bore Columns		105	10	
	TR-L05010	Cup Liner	453 00320	105	1	
	TR-L05010-5	Cup Liner		105	5	
	TR-L05010-10	Cup Liner		105	10	
	TR-L05012	PTV Liner, 2 mm ID	453 22045	120	1	
	TR-L05012-5	PTV Liner 2 mm ID		120	5	
	TR-L05012-10	PTV Liner 2 mm ID		120	10	

Liners for Shimadzu Capillary GCs


	Cat.No	Description	Similar Item No.	OD (mm)	Lenght (mm)	Pk
	TR-L05020	Split Liner 14A	221-32544-01	99	99	1
	TR-L05020-5	Split Liner 14A		99	99	5
	TR-L05020-10	Split Liner 14A		99	99	10
	TR-L05022	Splitless Liner 14A	221-32544-00	99	99	1
	TR-L05022-5	Splitless Liner 14A		99	99	5
	TR-L05022-10	Splitless Liner 14A		99	99	10

Liners for Varian Capillary GCs for Injector 1177



	Cat.No	Description	Similar Item No.	OD (mm)	Lenght (mm)	Pk
	TR-L05030	2 mm Splitless		6.5	79	1
	TR-L05030-5	2 mm Splitless		6.5	79	5
	TR-L05030-10	2 mm Splitless		6.5	79	10
	TR-L05032	2 mm Splitless, Quartz	392611924	6.5	79	1
	TR-L05032-5	2 mm Splitless, Quartz		6.5	79	5
	TR-L05032-10	2 mm Splitless, Quartz		6.5	79	10
	TR-L05034	4 mm Split/Splitless		6.3	79	1
	TR-L05034-5	4 mm Split/Splitless		6.3	79	5
	TR-L05034-10	4 mm Split/Splitless		6.3	79	10
	TR-L05036	4 mm Split/Splitless, Quartz		6.3	79	1
	TR-L05036-5	4 mm Split/Splitless, Quartz		6.3	79	5
	TR-L05036-10	4 mm Split/Splitless, Quartz		6.3	79	10
	TR-L05038	4 mm Split/Splitless, with Deactivated Glass Wool	392611934	6.3	79	1
	TR-L05038-5	4 mm Split/Splitless, with Deactivated Glass Wool		6.3	79	5
	TR-L05038-10	4 mm Split/Splitless, with Deactivated Glass Wool		6.3	79	10
	TR-L05040	Split Liner 4 mm id with Cup	392611931	6.3	79	1
	TR-L05040-5	Split Liner 4 mm id with Cup		6.3	79	5
	TR-L05040-10	Split Liner 4 mm id with Cup		6.3	79	10
	TR-L05042	Split Cup Liner, 4 mm id with Deactivated Glass Wool	392611932	6.3	79	1
	TR-L05042-5	Split Cup Liner, 4 mm id with Deactivated Glass Wool		6.3	79	5
	TR-L05042-10	Split Cup Liner, 4 mm id with Deactivated Glass Wool		6.3	79	10
	TR-L05044	Split Cup Liner, 4 mm id with OV-1/Chrom. W-P 80/100	392611933	6.3	79	1
	TR-L05044-5	Split Cup Liner, 4 mm id with OV-1/Chromosorb W-P 80/100		6.3	79	5
	TR-L05044-10	Split Cup Liner, 4 mm id with OV-1/Chromosorb W-P 80/100		6.3	79	10
	TR-L05046	Single Taper Liner, 2 mm id	392611926	6.5	79	1
	TR-L05046-5	Single Taper Liner, 2 mm id		6.5	79	5
	TR-L05046-10	Single Taper Liner, 2 mm id		6.5	79	10
	TR-L05048	Single Taper Liner, 4 mm id	392611927	6.5	79	1
	TR-L05048-5	Single Taper Liner, 4 mm id		6.5	79	5
	TR-L05048-10	Single Taper Liner, 4 mm id		6.5	79	10

	Cat.No	Description	Similar Item No.	OD (mm)	Lenght (mm)	Pk
	TR-L05050	Single Taper Liner, 4 mm id with Deactivated Glass Wool	392611936	6.5	79	1
	TR-L05050-5	Single Taper Liner, 4 mm id with Deactivated Glass Wool		6.5	79	5
	TR-L05050-10	Single Taper Liner, 4 mm id with Deactivated Glass Wool		6.5	79	10
	TR-L05052	Double Taper Liner, 4 mm id	392611929	6.5	79	1
	TR-L05052-5	Double Taper Liner, 4 mm id		6.5	79	5
	TR-L05052-10	Double Taper Liner, 4 mm id		6.5	79	10


Liners for Varian Capillary GCs for Injector 1075/1077

	Cat.No	Description	Similar Item No.	OD (mm)	Lenght (mm)	Pk
	TR-L05054	4 mm Open Split Liner	16-000830-00		72	1
	TR-L05054-5	4 mm Open Split Liner			72	5
	TR-L05054-10	4 mm Open Split Liner			72	10
	TR-L05056	4 mm Open Split Liner with Glass Wool	01-900109-01		72	1
	TR-L05056-5	4 mm Open Split Liner with Glass Wool			72	5
	TR-L05056-10	4 mm Open Split Liner with Glass Wool			72	10
	TR-L05058	Frit Split Liner			72	1
	TR-L05058-5	Frit Split Liner	01-900109-3		72	5
	TR-L05058-10	Frit Split Liner	16-000830-01		72	10
	TR-L05060	Splitless, Borosilicate Glass	01-900109-05		74	1
	TR-L05060-5	Splitless, Borosilicate Glass	03-949437-90		74	5
	TR-L05060-10	Splitless, Borosilicate Glass	03-949437-00		74	10
	TR-L05062	Splitless, Quartz			74	1
	TR-L05062-5	Splitless, Quartz			74	5
	TR-L05062-10	Splitless, Quartz			74	10

Liners for Varian Capillary GCs for Injector 1093/1094

	Cat.No	Description	Similar Item No.	OD (mm)	Lenght (mm)	Pk
	TR-L05064	0.5 mm SPI			54	1
	TR-L05064-5	0.5 mm SPI	01-900109-06		54	5
	TR-L05064-10	0.5 mm SPI	03-918332-01		54	10
	TR-L05066	0.8 mm SPI			54	1
	TR-L05066-5	0.8 mm SPI	01-900109-07		54	5
	TR-L05066-10	0.8 mm SPI	03-918332-02		54	10

Liners for Varian Capillary GCs for Injector 1078/1079

	Cat.No	Description	Similar Item No.	OD (mm)	Lenght (mm)	Pk
	TR-L05068	3.4 mm ID X 5 mm OD			54	1
	TR-L05068-5	3.4 mm ID X 5 mm OD			54	5
	TR-L05068-10	3.4 mm ID X 5 mm OD			54	10