

Low Volume Vessels, Conversion Kits



Agilent TruCenter 100ml suited for 7000, 7010, 7020 and 7025 Dissolution System.

Part No.	Description	OEM Ref
PSKIT100-VK	Kit includes;	12-0308
PSGLA100-01	100ml TruCenter Clear glass vessel	12-5145; 52-5010
PSCOVERV-LL200	Low evaporation cover with Cannula Hole	-
PSSMVASSY-TC	Small Volume Magnet. Adapter TruCenter Vessel	-
PSMINEPD019-01	Spin on/off mini paddle, 3165.S.	12-1443; 13-3608



PSKIT100-VK

Agilent TruCenter 200ml suited for 7000, 7010, 7020 and 7025 Dissolution System.

Part No.	Description	OEM Ref
PSKIT200-VK	Kit includes;	12-0312
PSGLA200-VK	200ml TruCenter Clear glass vessel	12-5033
PSCOVERV-LL200	Low evaporation cover with Cannula Hole	-
PSSMVASSY-TC	Small Volume Magnet. Adapter TruCenter Vessel	-
PSMINEPD019-01	Spin on/off mini paddle, 3165.S.	12-1443; 13-3608



PSKIT200-VK

Agilent Agilent EaseAlign 100ml suited for 705-DS, 7000 and 7010 Dissolution System.

Part No.	Description	OEM Ref
PSKIT100-01	Kit includes;	
PSGLA100-01	100ml EaseAlign Clear glass vessel	12-5040; 12-5145; 52-5010
PSSMVCOV-VK	Low evaporation cover with Cannula Hole	-
PSSMVASSY-VK	Small Volume Vessel Conversion Adapter, incl. cover	12-6050
PSMINEPD-01	Spin on/off mini paddle, 3165.S.	12-1443; 13-3608



PSKIT100-01

Agilent Agilent EaseAlign 200ml suited for 705-DS, 7000 and 7010 Dissolution System.

Part No.	Description	OEM Ref
PSKIT200-01	Kit includes;	
PSGLA200-01	200ml EaseAlign Clear glass vessel	12-5050
PSSMVCOV-VK	Low evaporation cover with Cannula Hole	-
PSSMVASSY-VK	Small Volume Vessel Conversion Adapter, incl. cover	12-6050
PSMINEPD-01	Spin on/off mini paddle, 3165.S.	12-1443; 13-3608



PSKIT200-01

Low Volume Vessels, Conversion Kits

Agilent TruAlign 100ml suited for 708-DS Dissolution System

Part No.	Description	OEM Ref
PSKIT100-TA	Kit includes;	
PSGLA100-TA	100ml TruAlign Clear glass vessel	12-5145
PSCOVERV-LL200	Low evaporation cover with Cannula Hole	-
PSSMVASSY-TA	Small Volume Vessel Conversion Adapter	-
PSMINEPD-01	Spin on/off mini paddle, 316S.S.	12-1443; 13-3608



PSKIT100-TA

Agilent TruAlign 200ml suited for 708-DS Dissolution System.

Part No.	Description	OEM Ref
PSKIT200-TA	Kit includes;	
PSGLA200-TA	200ml Clear glass vessel	12-5147
PSCOVERV-LL200	Low evaporation cover with Cannula Hole	-
PSSMVASSY-TA	Small Volume Vessel Conversion Adapter	-
PSMINEPD-01	Spin on/off mini paddle, 316S.S.	12-1443; 13-3608



PSKIT200-TA

Distek 100ml suited for 2100, 2500 and Evolution 6300 Dissolution System.

Part No.	Description	OEM Ref
PSKIT100-DK	Kit includes;	2800-0600
PSGLA100-DK	100ml Clear glass vessel	3010-0094
PSSMVCOV-DK	Low evaporation cover with Cannula Hole	-
PSSMVASSY-DK	Small Volume Vessel Conversion Adapter, incl. cover	2800-0602
PSMINBLD-DK	Detachable mini paddle blade, 316S.S.	3200-0185



PSKIT100-DK

Distek 200ml suited for 2100, 2500 and Evolution 6300 Dissolution System.

Part No.	Description	OEM Ref
PSKIT200-DK	Kit includes;	
PSGLA200-DK	200ml Clear glass vessel	3100-0095
PSSMVCOV-DK	Low evaporation cover with Cannula Hole	-
PSSMVASSY-DK	Small Volume Vessel Conversion Adapter, incl. cover	2800-0602
PSMINBLD-DK	Detachable mini paddle blade, 316S.S.	3200-0185



PSKIT200-TA



Erweka 100ml LH suited for DT126, DT720, DT820, DT1410 and DT1610 Dissolution System.

Part No.	Description	OEM Ref
PSKIT100-EWLH	Kit includes;	
PSGLA100-EW	100ml Clear glass vessel	20482
PSSMVASSY-EW	Small Volume Vessel Conversion Adapter	-
PSCOVERV-LL200	Low evaporation cover with Cannula Hole	-
PSMINEPD-EWLH	Spin on/off mini paddle, 316S.S., LH	-

Erweka 100ml HH suited for DT126, DT720, DT820, DT1410 and DT1610 Dissolution System.

Part No.	Description	OEM Ref
PSKIT100-EWHH	Kit includes;	
PSGLA100-EW	100ml Clear glass vessel	20482
PSSMVASSY-EW	Small Volume Vessel Conversion Adapter	-
PSCOVERV-LL200	Low evaporation cover with Cannula Hole	-
PSMINEPD-EWHH	Spin on/off mini paddle, 316S.S., HH	-

Erweka 400ml LH suited for DT126, DT720, DT820, DT1410 and DT1610 Dissolution System.

Part No.	Description	OEM Ref
PSKIT400-EWLH	Kit includes;	22399
PSGLA400-EW	400ml Clear glass vessel	18373
PSSMVASSY-EW	Small Volume Vessel Conversion Adapter	-
PSCOVERV-LL400	Low evaporation cover with Cannula Hole	-
PSMINEPD-EWLH	Spin on/off mini paddle, 316S.S., LH	-

Erweka 400ml HH suited for DT126, DT720, DT820, DT1410 and DT1610 Dissolution System.

Part No.	Description	OEM Ref
PSKIT400-EWHH	Kit includes;	22398
PSGLA400-EW	400ml Clear glass vessel	18373
PSSMVASSY-EW	Small Volume Vessel Conversion Adapter	-
PSCOVERV-LL400	Low evaporation cover with Cannula Hole	-
PSMINEPD-EWHH	Spin on/off mini paddle, 316S.S., HH	-



Low Volume Vessels, Conversion Kits

Hanson 150ml suited for SR6, SR8-Plus Dissolution System.

Part No.	Description	OEM Ref
PSKIT150-HR	Kit includes;	
PSGLA150R-HR	150ml Clear glass vessel, round bottom	74-105-151; 72-600-573
PSCOVERH-LL200	Low evaporation cover with Cannula Hole	72-600-560
PSSMVASSY-HR	Small Volume Vessel Conversion Adapter	-
PSMINEPD-HR	Spin on/off mini paddle, 316S.S.	74-105-203; 74-105-204



PSKIT150-HR



PSMINEPD-HR

Hanson 150ml suited for Vision Dissolution System.

Part No.	Description	OEM Ref
PSKIT150-HRV	Kit includes;	
PSGLA150R-HRV	150ml Clear glass vessel, round bottom	74-105-151; 72-600-573
PSCOVERH-LL200	Low evaporation cover with Cannula Hole	72-600-560
PSSMVASSY-HRV	Small Volume Vessel Easi-Lock Adapter	74-105-155
PSMINEPD-HR	Spin on/off mini paddle, 316S.S.	74-105-203; 74-105-204



PSKIT150-HRV



PSMINEPD-HR

Sotax Mini vessels for AT systems.

Part No.	Description	OEM Ref
PSGLA200-ST	Clear glass mini vessel, round bottom	14535-01
PSGLA200-AST	Amber glass mini vessel, round bottom	-
PSGLA200F-ST	Clear glass mini vessel, flat bottom	14534-01