

# Standard Vessels



## Agilent TruCenter Compatible Vessels

Part No.	Description	OEM Ref
PSGLA100-VK	100 ml TruCenter, Clear Glass Vessel	12-5031; 12-5040
PSGLA100-AVK	100 ml TruCenter, Amber Glass Vessel	-
PSGLA200-VK	200 ml TruCenter, Clear Glass Vessel	12-5033
PSGLA200-AVK	200 ml TruCenter, Amber Glass Vessel	-
PSGLA900-VK	1000 ml TruCenter, Clear Glass Vessel for 7025 and 7030	-
PSGLA900-AVK	1000 ml TruCenter, Amber Glass Vessel for 7025 and 7030	12-5020
PSGLA900-VKC	1000 ml TruCenter, Clear Glass Vessel with Collar for 7025 and 7030	12-5035
PSGLA900-AVKC	1000 ml TruCenter, Amber Glass Vessel with Collar for 7025 and 7030	12-5120
PSGLA9PK-VK	1000 ml TruCenter Clear Glass Vessel, Peak, for V-Series Bath	12-5125
PSGLA9PK-AVK	1000 ml TruCenter Amber Glass Vessel, Peak, for V-Series Bath	12-5505
PSGLA9PK-VKC	1000 ml TruCenter Clear Glass Vessel, Peak, with Collar, for V-Series Bath	12-5125
PSGLA9PK-AVKC	1000 ml TruCenter Amber Glass Vessel, Peak, with Collar, for V-Series Bath	12-5130
PSGLA02K-VKC	2000 ml TruCenter, Clear Glass Vessel, with Collar	-
PSGLA02K-AVKC	2000 ml TruCenter, Amber Glass Vessel, with Collar	-
PSGLA2KPK-VKC	2000 ml TruCenter, Clear Glass Vessel, Peak with Collar	-
PSGLA2KPK-AVKC	2000 ml TruCenter, Amber Glass Vessel, Peak with Collar	-

Collar Centering Ring for TruCenter Vessel  
PSCNTRNG-VK



PSGLA900-VKC



PSGLA900-VK



PSGLA900-AVK



PSGLA9PK-VK



PSGLA9PK-AVK



PSGLA02K-VKC



PSGLA02K-AVKC

## Agilent EaseAlign Compatible Vessels

Part No.	Description	OEM Ref
PSGLA100-01	100 ml, EaseAlign, Clear Glass Vessel	12-5040
PSGLA100-A01	100 ml, EaseAlign, Amber Glass Vessel	12-5041
PSGLA200-01	200 ml, EaseAlign, Clear Glass Vessel	12-5050
PSGLA200-A01	200 ml, EaseAlign, Amber Glass Vessel	12-5051
PSGLA500-01	500 ml, EaseAlign, Clear Glass Vessel	12-5060
PSGLA500-A01	500 ml, EaseAlign, Amber Glass Vessel	12-5065
PSGLA900-01	1000 ml, EaseAlign, Clear Glass Vessel	12-5000
PSGLA900-A01	1000 ml, EaseAlign, Amber Glass Vessel	12-5010
PSGLA900-01T	1000 ml, EaseAlign, Clear Glass Vessel, Teflon Coated	12-5020
PSGLA9PK-01	1000 ml, EaseAlign, Clear Glass Vessel, Peak	12-5500
PSGLA9PK-A01	1000 ml, EaseAlign, Amber Glass Vessel, Peak	12-5505
PSGLA9PK-01T	1000 ml, EaseAlign, Clear Glass Vessel, Peak, Teflon coated	12-5525
PSPLA900-01	1000 ml Clear Plastic Vessel	12-5200
PSPLA900-A01	1000 ml Amber Plastic Vessel	12-5210
PSGLA02K-01	2000 ml, EaseAlign, Clear Glass Vessel	12-5070
PSGLA02K-A01	2000 ml, EaseAlign, Amber Glass Vessel	12-5075
PSGLA2KPK-01	2000 ml, EaseAlign, Clear Glass Vessel, Peak	12-5510
PSGLA2KPK-A01	2000 ml, EaseAlign, Amber Glass Vessel, Peak	-

# Standard Vessels

## Agilent EaseAlign Compatible Vessels

Part No.	Description	OEM Ref
PSEZLBKT-01	Vessel Centering Bracket, 6 Replaceable Pins	12-6050
PSEZLKIT-01	Vessel Centering Bracket Repair Kit (Set of 12 pins)	12-6065
PSEZLPEG-01	Vessel Centering Bracket Peg	12-6065

Centering bracket



PSEZLBKT-01

## Agilent TruAlign Compatible Vessels

Part No.	Description	OEM Ref
PSGLA100-TA	100 ml TruAlign, Clear Glass Vessel	12-5145
PSGLA100-ATA	100 ml TruAlign, Amber Glass Vessel	12-5146
PSGLA200F-TA	200 ml TruAlign, Clear Glass Vessel, Flat Bottom	12-5170
PSGLA200F-ATA	200 ml TruAlign, Amber Glass Vessel, Flat Bottom	-
PSGLA200-TA	200 ml TruAlign, Clear Glass Vessel	12-5147
PSGLA200-ATA	200 ml TruAlign, Amber Glass Vessel	-
PSCSV250-TA	250 ml TruAlign, Clear Glass Vessel, Chinese Pharmacopeia	12-1452
PSCSV250-ATA	250 ml TruAlign, Amber Glass Vessel, Chinese Pharmacopeia	-
PSGLA900-TA	1000 ml TruAlign, Clear Glass Vessel	12-5149
PSGLA900-ATA	1000 ml TruAlign, Amber Glass Vessel	12-5152
PSGLA9PK-TA	1000 ml TruAlign, Clear Glass Vessel, Peak	12-5153
PSGLA9PK-ATA	1000 ml TruAlign, Amber Glass Vessel, Peak	-
PSGLA02K-TA	2000 ml TruAlign, Clear Glass Vessel	12-5157
PSGLA02K-ATA	2000 ml TruAlign, Amber Glass Vessel	12-5158
PSGLA2KPK-TA	2000 ml TruAlign, Clear Glass Vessel, Peak	12-5129
PSGLA2KPK-ATA	2000 ml TruAlign, Amber Glass Vessel, Peak	-

TruAlign centering ring mounted to the vessel

PSGLA900-TA



## Caleva Compatible Vessels

Part No.	Description	OEM Ref
PSGLA900-CA	1000ml Clear Glass Vessel	0810004
PSGLA900-ACA	1000ml Amber Glass Vessel	0810009
PSGLA9PK-CA	1000 ml Clear Glass Vessel, Peak	-
PSGLA9PK-ACA	1000 ml Amber Glass Vessel, Peak	-
PSGLA02K-CA	2000 ml Clear Glass Vessel	0810005
PSGLA02K-ACA	2000 ml Amber Glass Vessel	0810010
PSGLA2KPK-CA	2000 ml Clear Glass Vessel, Peak	-
PSGLA2KPK-ACA	2000 ml Amber Glass Vessel, Peak	-



PSGLA900-CA



PSGLA900-ACA

## Copley Compatible Vessels

Part No.	Description	OEM Ref
PSGLA900-CP	1000ml Clear Glass Vessel	-
PSGLA900-ACP	1000ml Amber Glass Vessel	-
PSGLA9PK-CP	1000 ml Clear Glass Vessel, Peak	-
PSGLA9PK-ACP	1000 ml Amber Glass Vessel, Peak	-



PSGLA900-CO



PSGLA900-ACO



## Distek Compatible Vessels

Part No.	Description	OEM Ref
PSGLA200-DK	200 ml Clear Glass Vessel	3100-0095
PSGLA200-ADK	200 ml Amber Glass Vessel	3100-0095-AMB
PSCSV250-DK	250 ml Clear Glass Vessel, Chinese Pharmacopeia	3010-0090
PSCSV250-ADK	250 ml Amber Glass Vessel, Chinese Pharmacopeia	–
PSGLA500-DK	500 ml Clear Glass Vessel	3010-0040
PSGLA500-ADK	500 ml Amber Glass Vessel	3010-0040-AMB
PSGLA900-DK	1000 ml Clear Glass Vessel	3010-0069
PSGLA900-ADK	1000 ml Amber Glass Vessel	3010-0023
PSGLA900F-DK	1000ml Clear Glass Vessel, Flat bottom	3010-0022
PSGLA905F-ADK	1000ml Amber Glass Vessel, Flat bottom	–
PSGLA9CR-DK	1000 ml Clear Glass Vessel with Acculign Ring	3010-0093
PSGLA9CR-ADK	1000 ml Amber Glass Vessel with Acculign ring	3010-0093-AMB
PSGLA9PK-DK	1000 ml Clear Glass Vessel, Peak	3010-0550
PSGLA9PK-ADK	1000 ml Amber Glass Vessel, Peak	3010-0550-AMB
PSGLA900-DKT	1000 ml Clear Glass Vessel, Teflon Coated with Acculign Ring	3010-0002-T
PSGLP900-DK	1000 ml Clear Glass Vessel, Premier	3010-0096
PSGLP900-ADK	1000 ml Amber Glass Vessel, Premier	3010-0096-AMB
PSGLP900-DKT	1000 ml Clear Glass Vessel, Premier, Teflon Coated	3010-0096-T
PSPLA900-DK	1000 ml Clear Plastic Vessel	3350-0095
PSPLA900-ADK	1000 ml Amber Plastic Vessel	3350-0095-AMB
PSGLA02K-DK	2000 ml Clear Glass Vessel with Acculign Ring	3010-0285
PSGLA02K-ADK	2000 ml Amber Glass Vessel with Acculign Ring	3010-0327-AMB
PSGLA2KPK-DK	2000 ml Clear Glass Vessel, Peak	–
PSGLA2KPK-ADK	2000 ml Amber Glass Vessel, Peak	–
PSGLA04K-DK	4000 ml Clear Glass Vessel with Acculign Ring	3010-0286
PSGLA04K-ADK	4000 ml Amber Glass Vessel with Acculign Ring	3010-0286-AMB
PSGLA4KPK-DK	4000 ml Clear Glass Vessel, Peak	–
PSGLA4KPK-ADK	4000 ml Amber Glass Vessel, Peak	–
PSRING01-DK	(Acculign Ring) Plastic Centering Ring	3350-0096
PSVSCCLIP-DK	Vessel Hold Down Clip (for 6300 and 2100 Series)	2821-0257



PSGLA900-DK



PSGLA900-ADK



PSGLA900F-DK



PSGLA905F-ADK



PSGLA9PK-DK



PSGLA9CR-ADK



PSGLA02K-DK



PSGLA500-DK



PSGLA200-DK



PSRING01-DK



PSVSCCLIP-DK

# Standard Vessels

## Electrolab Compatible Vessels

Part No.	Description	OEM Ref
PSGLA100-EL	100ml Clear Glass Vessel	0103A00001
PSGLA250-EL	250ml Clear Glass Vessel	0103A00002
PSGLA900-EL	1000 ml Clear Glass Vessel	0103A00004
PSGLA900-AEL	1000 ml Amber Glass Vessel	0103A00005
PSGLA9PK-EL	1000 ml Clear Glass Vessel, Peak	0103A00010
PSGLA9PK-AEL	1000 ml Amber Glass Vessel, Peak	0103A00018
PSPLA900-EL	1000ml Clear Plastic Vessel	0103A00038
PSPLA900-AEL	1000ml Amber Plastic Vessel	0103A00039
PSGLA02K-EL	2000 ml Clear Glass Vessel	0103A00006
PSGLA02K-AEL	2000 ml Amber Glass Vessel	-
PSGLA2KPK-EL	2000 ml Clear Glass Vessel, Peak	-
PSGLA2KPK-AEL	2000 ml Amber Glass Vessel, Peak	-



PSGLA900-EL



PSGLA900-AEL



PSGLA900-EW



PSGLA900-AEW

## Erweka Compatible Vessels

Part No.	Description	OEM Ref
PSGLA250-EW	250 ml Clear Glass Vessel	80-000-0202
PSGLA250-AEW	250 ml Amber Glass Vessel	-
PSCVS250-EW	250 ml Clear Glass Vessel, Chinese Pharmacopeia	-
PSCVS250-AEW	250 ml Amber Glass Vessel, Chinese Pharmacopeia	-
PSGLA900-EW	1000 ml Clear Glass Vessel	80-000-1002; 18365
PSGLA900-AEW	1000 ml Amber Glass Vessel	80-000-1004; 18366
PSGLA9PK-EW	1000 ml Clear Glass Vessel, Peak	19115
PSGLA9PK-AEW	1000 ml Amber Glass Vessel, Peak	-
PSPLA900-EW	1000 ml Clear Plastic Vessel	80-000-1001
PSPLA900-AEW	1000 ml Amber Plastic Vessel	-
PSGLA02K-EW	2000 ml Clear Glass Vessel	80-000-2002; 18367
PSGLA02K-AEW	2000 ml Amber Glass Vessel	18368
PSCNTRNG-EW	Centering Ring for Erweka Vessel	84-060-3000



PSGLA02K-EW



PSPLA900-EW



PSCNTRNG-EW



## Hanson Compatible Vessels SR6-SR8

Part No.	Description	OEM Ref
PSGLA150R-HR	150 ml Clear Glass Vessel	72-600-573
PSGLA150R-AHR	150 ml Amber Glass Vessel	72-600-583
PSGLA150F-HR	150 ml Clear Glass Vessel, Flat Bottom	74-105-152
PSGLA150F-AHR	150 ml Amber Glass Vessel, Flat Bottom	74-105-162
PSCSV250-HR	250 ml Clear Glass, Chinese Pharmacopeia	74-105-280
PSCSV250-AHR	250 ml Amber Glass, Chinese Pharmacopeia	74-105-281
PSGLA500-HR	500 ml Clear Glass Vessel	-
PSGLA500-AHR	500 ml Amber Glass Vessel	-
PSGLA900-HR	1000 ml Clear Glass Vessel, USP, Precision	72-600-510
PSGLA900-AHR	1000 ml Amber Glass Vessel, USP, Precision	72-600-556
PSPLA900-HR	1000 ml Plastic Clear Vessel, USP, Precision	91-030-003
PSPLA900-AHR	1000 ml Plastic Amber Vessel, USP, Precision,	91-030-004
PSGLA9PK-HR	1000 ml Clear Glass Vessel, Peak, USP, Precision	95-590-047
PSGLA9PK-AHR	1000 ml Amber Glass Vessel, Peak, USP, Precision	95-590-048
PSGLA02K-HR	2000 ml Clear Glass Vessel, USP, Precision	72-600-509
PSGLA02K-AHR	2000 ml Amber Glass Vessel, USP, Precision	72-600-580
PSGLA2KPK-HR	2000 ml Clear Glass Vessel, Peak	-
PSGLA2KPK-AHR	2000 ml Amber Glass Vessel, Peak	-
PSGLA04K-HR	4000 ml Clear Glass Vessel	-
PSCNTRNG-HR	Vessel Centering Ring (1 & 2 L Vessels)	73-100-550



## Hanson Compatible Vessels for Vision

Part No.	Description	OEM Ref
PSGLA900-HRV	1000ml Clear Glass Vessel, Vision Easi-lock Comp., USP Precision	74-104-101
PSGLA900-AHRV	1000ml Amber Glass Vessel, Vision Easi-lock Comp., USP Precision	74-104-102
PSGLA9PK-HRV	1000ml Clear Glass PEAK Vessel, Vision Easi-lock Comp., USP Precision	74-104-131
PSGLA9PK-AHRV	1000ml Amber Glass PEAK Vessel, Vision Easi-lock Comp., USP Precision	74-104-132
PSGLA02K-HRV	2000ml Clear Glass Vessel, Vision Easi-lock Comp., USP Precision	74-104-111
PSGLA02K-AHRV	2000ml Amber Glass Vessel, Vision Easi-lock Comp., USP Precision	74-104-112
PSCNTRNG-HRV	Vessel centering Ring (1L & 2L Vessels), Vision	74-105-106



# Standard Vessels

## Logan Compatible Vessels

Part No.	Description	OEM Ref
PSGLA900-LG	1000ml Clear Glass Vessel	200-2045; 800-0124
PSGLA9PK-LG	1000ml Clear Glass Vessel, PEAK	800-0140
PSGLA900-ALG	1000ml Amber Glass Vessel	800-0147
PSPLA900-LG	1000ml Clear Plastic Vessel	200-0900
PSPLA900-ALG	1000ml Amber Plastic Vessel	–
PSGLA02K-LG	2000ml Clear Glass Vessel	200-0126
PSGLA02K-ALG	2000ml Amber Glass Vessel	–



PSGLA900-LG



PSGLA900-ALG

## Pharmatest Compatible Vessels

Part No.	Description	OEM Ref
PSGLA900-PT	1000 ml Clear Glass Vessel	303-5225
PSGLA900-APT	1000 ml Amber Glass Vessel	303-5425
PSGLA900-PTR	1000 ml Clear Glass Vessel, Enlarged Top Flange	–
PSGLA900-APTR	1000 ml Amber Glass Vessel, Enlarged Top Flange	–
PSGLA9PK-PT	1000 ml Clear Glass Vessel, Peak	–
PSGLA9PK-APT	1000 ml Amber Glass Vessel, Peak	–
PSPLA900-PT	1000 ml Clear Plastic Vessel	–
PSPLA900-APT	1000 ml Amber Plastic Vessel	–
PSGLA02K-PT	2000 ml Clear Glass Vessel	303-5221
PSGLA02K-APT	2000 ml Amber Glass Vessel	303-5222
PSGLA2KPK-PT	2000 ml Clear Glass Vessel, Peak	–
PSGLA2KPK-APT	2000 ml Amber Glass Vessel, Peak	–



PSGLA900-PT



PSGLA900-APT



PSPLA900-PT



PSPLA900-APT

## Sotax Compatible Vessels

Part No.	Description	OEM Ref
PSCSV250-ST	250 ml mini Chinese Clear Glass Vessel	9650
PSCSV250-AST	250 ml mini Chinese Amber Glass Vessel	9650-2
PSGLA900-ST	1000 ml Clear Glass Vessel	2503
PSGLA900-AST	1000 ml Amber Glass Vessel	2503-2



PSGLA900-ST



PSCSV250-ST

## Sotax Compatible Vessels Xtend

Part No.	Description	OEM Ref
PSCSV250-STX	250ml Clear Glass Chinese Small Volume Vessel, Xtend	14536-01
PSCSV250-ASTX	250ml Amber Glass Chinese Small Volume Vessel, Xtend	–
PSGLA900-STX	1000ml Clear Glass Vessel, Xtend	14143-01
PSGLA900-ASTX	1000ml Amber Glass Vessel, Xtend	14143-02



PSCSV250-STX



PSGLA900-ASTX



## Toyama Compatible Vessels

Part No.	Description	OEM Ref
PSGLA900-TO	1000ml Clear Glass Vessel	-
PSGLA900-ATO	1000ml Amber Glass Vessel	-
PSGLA9PK-TO	1000 ml Clear Glass Vessel, Peak	-
PSGLA9PK-ATO	1000 ml Amber Glass Vessel, Peak	-



PSGLA900-TO



PSGLA900-ATO

## Zymark Compatible Vessels

Part No.	Description	OEM Ref
PSGLA900-ZM	1000 ml Clear Glass Vessel	-
PSGLA900-AZM	1000 ml Amber Glass Vessel	-
PSGLA9CR-ZM	1000 ml Clear Glass Vessel with Centering Ring	58730
PSGLA9CR-AZM	1000 ml Amber Glass Vessel with Centering Ring	-
PSGLA9PK-ZM	1000 ml Clear Glass Vessel, Peak	-
PSGLA9PK-AZM	1000 ml Amber Glass Vessel, Peak	-
PSCNTRNG-ZM	Multidose Centering Ring	69606



PSGLA900-ZM



PSGLA900-AZM

## Flange Protection Ring

Part No.	Description	OEM Ref
PSPRORNG-01	Protection ring for Agilent EaseAlign 1000ml vessel	-
PSPRORNG-DK	Protection ring for Distek 1000ml vessel	-
PSPRORNG-HR	Protection ring for Hanson SR6/SR8 1000ml vessel	-
PSPRORNG-PT	Protection ring for Pharmatest 1000ml vessel	-
PSPRORNG-ST	Protection ring for Sotax AT7 1000ml vessel	-



PSPRORNG-01

### Teflon Coated Vessels are available on customers request.

In many cases Teflon Coated Vessels are used when the measured sample is sticky and sticks to the wall or bottom of the dissolution vessel. This can affect your measurement dramatically. The Teflon coating on the inside of the dissolution vessel can prevent that the sample sticks to the dissolution vessel.




## Universal Vessel Stands

For the safe transportation and storage of vessels

Due to their shape, dissolution vessels are difficult to hold securely. This can often result in scratched, chipped or even smashed vessels.

The steel/acrylic stands are available in 6 or 8 vessel capacities and can hold vessels from any dissolution bath manufacturer.

They can also be stacked together to save bench space and enable you to transport many vessels at a time.



Part No.	Description	OEM Ref
PSVESSTD06-UN	Vessel Stand, 6 Position	-
PSVESSTD08-UN	Vessel Stand, 8 Position	-
PSVESSTD08-ST	Vessel Stand Sotax, 8 Positions	-