

# Misc Accessories

## Tubing

| Part No.    | Description  | OEM Ref   |          |
|-------------|--|-----------|----------|
| PSPPTUBV-01 | Peristaltic Pump Tubing, Santoprene, Coded Purple & White (Set of 8)                             | 17-6000   | <b>A</b> |
| PSPPTUBV-02 | Peristaltic Pump Tubing, Silicone, Coded Purple & White (Set of 8)                               | 17-6010   | <b>A</b> |
| PSPPTUBD-01 | Peristaltic Pump Tubing, Santoprene (Set of 6)   | 5640-0337 | <b>D</b> |
| PSTAATUB-HP | 3 Collar Tygon Auto Analyser Tubing, 0.08" (2.06mm) ID (Set of 12)                               | 5041-2166 | <b>D</b> |
| PSTEFTUB-HP | FEP PTFE Tubing for HP 1/32" (0.8mm) ID x 1/16" (1.6mm) OD, 500' (152m) Roll                     | -         |          |
| PSTEFTUB-8V | Colour Coded Sample PTFE Tubing, 0.42" (10.67mm) ID x 0.066" (1.68mm) OD, 8 Lengths of 50' (15m) | 17-6035   | <b>A</b> |
| PSTUBING    | Tygon Tubing, 1/16" (1.6mm) x 1/8" (3.2mm) (15m)   | -         |          |



PSTEFTUB-8V



PSPPTUBV-01



PSPPTUBV-02



PSTAATUB-HP



PSTUBING



PSTEFTUB-HP



PSPPTUBD-02

**D** = Distek    **A** = Agilent

## Fittings, Tubing & Couplings

| Part No.     | Description   | OEM Ref    |
|--------------|---|------------|
| PS91-601-130 | Fitting, 1/16 in. Ferrule, Blue , pk/10                             | -          |
| PS91-601-131 | Fitting, 1/4-28 x 1/16 in. Nut, Black, pk/10                        | -          |
| PS91-601-132 | Fitting, 1/4-28 x 1/16 in. Nut, Red, pk/10                          | -          |
| PS91-601-133 | Fitting, 1/4-28 x 1/16 in. Nut, White, pk/10                        | -          |
| PS91-601-134 | Fitting, 1/4-28 x 1/16 in. Nut, Yellow, pk/10                       | -          |
| PS91-601-135 | Fitting, 1/4-28 x 1/16 in. Nut, Green, pk/10                        | -          |
| PS91-601-136 | Fitting, 1/4-28 x 1/16 in. Nut, Blue, pk/10                         | -          |
| PS91-601-137 | Fitting, 1/4-28 x 1/16 in. Nut, Clear, pk/10                        | -          |
| PS91-601-138 | Fitting, 1/4-28 Union, White, pk/10                                 | -          |
| PS91-601-139 | Fitting, 1/4-28 x 1/16 in. Nut, Purple, pk/10                       | -          |
| PS91-601-140 | Fitting, 1/8 in. Ferrule, Yellow, pk/10                             | -          |
| PS91-601-141 | Fitting, 1/4-28 x 1/8 in. Nut, Black, pk/10                         | -          |
| PS91-601-142 | Fitting, 1/4-28 x 1/8 in. Nut, Red, pk/10                           | -          |
| PS91-601-143 | Fitting, 1/4-28 x 1/8 in. Nut, White, pk/10                         | -          |
| PS91-601-144 | Fitting, 1/4-28 x 1/8 in. Nut, Yellow, pk/10                        | -          |
| PS91-601-145 | Fitting, 1/4-28 x 1/8 in. Nut, Green, pk/10                         | -          |
| PS91-601-146 | Fitting, 1/4-28 x 1/8 in. Nut, Blue, pk/10                          | -          |
| PS91-601-147 | Fitting, 1/4-28 x 1/8 in. Nut, Clear, pk/10                         | -          |
| PS91-601-149 | Fitting, 1/4-28 x 1/8 in. Nut, Purple, pk/10                        | -          |
| PS91-601-261 | Fitting, 1/4-28 to Female Luer, MicroettePlus Cell Connector, pk/10 | -          |
| PS91-600-106 | Tubing, PTFE, 1/32 in. i.d. x 1/16 in. o.d. (100 ft. roll), pk/10   | -          |
| PS91-601-036 | Fitting, 1/4-28 Union, Clear, pk/10                                 | -          |
| PS91-601-172 | Fitting, 1/8-B x Female Luer, White, pk/10                          | 91-601-172 |
| PS91-601-176 | Fitting, 1/8-B x Male Luer, White, pk/10                            | 91-601-176 |
| PS91-601-237 | Fitting, 1/4" Malex1/4" Female, Stem Elbow, MMPlus, 10/pk           | 91-601-237 |
| PS91-601-262 | Fitting, female luer to female luer, pk/10                          | 91-601-262 |



Ferrule, blue, 1/16 in



PS91-601-131 PS91-601-132



PS91-601-133 PS91-601-134



PS91-601-135 PS91-601-136



PS91-601-137 PS91-601-139



91-601-261



PS91-601-237



PS91-601-262