

# Paddles

## Agilent Compatible Paddles

Part No.	Description	OEM Ref
PSEPD013-01	14.5" (370mm) Electropolished Paddle 316S.S. for 7025/7030 series	12-1035; 13-3590
PSEPD015-01	15" (380mm) Electropolished Paddle 316S.S. for 7000/7010 series	12-1030
PSEPD019-01	19" (485mm) Electropolished Paddle 316S.S. for 7000/7010/705/708/709-DS series	12-1040
PSEPD021-01	21" (530mm) Electropolished Paddle 316S.S. for 7000/7010/705/708/709-DS series	–
PSEPD024-01	24" (610mm) Electropolished Paddle 316S.S. for 7000/7010/705/708/709-DS series	12-1318; 12-1045; 13-3597
PSNPD013-01	14.5" (370mm) PTFE Coated Paddle 316S.S. for 7025/7030 series	12-1005; 13-3589
PSNPD015-01	15" (380mm) PTFE Coated Paddle 316S.S. for 7000/7010 series	12-1000
PSNPD019-01	19" (485mm) PTFE Coated Paddle 316S.S. or 7000/7010/705/708/709-DS series	12-1010



PSSHFT165-DK

## Distek Compatible Paddles

Part No.	Description	OEM Ref
PSEPD016-DK	16.5" (420mm) Electropolished Paddle, 316S.S.	2800-0239; 2800-0500
PSNPD016-DK	16.5" (420mm) PTFE Coated Paddle, 316S.S.	2821-0171
PSTPD016-DK	16.5" (420mm) Solid PTFE Paddle	2800-0178
PSPDLBLD-DK	Paddle Blade, screw on, 316S.S.	3251-1050
PSSHFT165-DK	16.5" (420mm) Detachable Shaft, 2100 series	2800-0165
PSTUMSCR-DK	Thumb Screw for Distek Compatible Paddle Shafts	2821-0094



PSPDLBLD-DK



PSTUMSCR-DK

## Electrolab Compatible Paddles

Part No.	Description	OEM Ref
PSEPD016-EL	Electropolished Paddle 316S.S.	0103A00085
PSEPD024-EL	Electropolished Paddle 316S.S. for 2L vessel	0103A00086
PSCSVEPD-EL	Electropolished Mini paddle 316S.S. for 250ml	0103A00084
PSMINEPD-EL	Electropolished Mini paddle 316S.S. for 100ml	0103A00083



PSSPNEPD-PT



PSPADHUB-EW

## Erweka Compatible Paddles

Part No.	Description	OEM Ref
PSEPD016-EW	16" (405mm) Electropolished Paddle 316S.S.	–
PSPADHUB-EW	Paddle blade Electropolished 316S.S.	81-000-7001; 22403



## Hanson Compatible Paddles

Part No.	Description	OEM Ref
PSEPD016-HR	16" (405mm) Electropolished Paddle, 316S.S.	65-700-036
PSNPD016-HR	16" (405mm) PTFE Coated Paddle, 316S.S.	65-700-030
PSTPD016-HR	16" (405mm) Solid PTFE Paddle	72-800-456
PSTPD018-HR	18" (460mm) Solid PTFE Paddle	72-800-461



PTFE Coated Paddle  
(PSNPD series)

## Logan Compatible Paddles

Part No.	Description	OEM Ref
PSSHFT175-LG	Detachable Shaft 17.5", 316S.S.	802-0043
PSSHFT023-LG	Detachable Shaft 23", 316S.S.	802-0101
PSPDLBLD-LG	Paddle Blade, screw on, 316S.S.	800-0059



Solid PTFE Paddle  
(PSTPD series)

## Pharmatest Compatible Paddles

Part No.	Description	OEM Ref
PSEPD610-PT	Electropolished Paddle 316S.S. for PTWS 610/1210	316-0210
PSSPNSHT14-PT	14" (355mm) Detachable Shaft, 316S.S.	316-0116
PSSPNEPD-PT	Paddle Adapter, 316S.S.	316-0210
PSSPNNPD-PT	Paddle Adapter, PTFE Coated	316-0250
PSCSVEPD-PT	Paddle Adapter for Chinese Small Volume	316SS




Polished Stainless Steel Paddle  
(PSEPD series)

## Sotax Compatible Paddles

Part No.	Description	OEM Ref
PSCSVEPD-ST	Electropolished Chinese Small volume Paddle 316S.S.	9656
PSEDPAD-AT7	Electropolished Paddle 316S.S.	7525
PSNPDAD-AT7	PTFE coated paddle	7525-2
PSPADHUB-ST	Paddle Electropolished 316 S.S with sample opening	2985-2
PSFILHUB-ST	Sieve for sampling opening, pk/10	2971



PSPADHUB-ST




PSEPDCSV-AT7

## Basket Shaft/Paddle Holders

Ideal for both basket shaft and paddle storage

### Key Benefits:

- Prevents shafts/paddles from becoming bent or scratched
- Allows quick and easy transport of many shafts/paddles
- Handy end plate to stop shafts/paddles sliding out when carrying
- Holds all lengths and types of basket shafts and paddles
- Non-corrosive, quality plastic construction
- Space-saving, attractive design
- Available in both 12 and 16 shaft capacities



**Dual Purpose:** This basket shaft/paddle holder is designed for optimal handling and storage of both basket shafts and paddles. Why not use one to keep a set of both paddles and basket shafts next to each bath?

Part No.	Description	OEM Ref
PSPADHLD-12A	Basket Shaft/Paddle Rack, 12 Position	-
PSPADHLD-16A	Basket Shaft/Paddle Rack, 16 Position	-
PSPADHLD-ST	Basket Shaft/Paddle Rack for Sotax, 16 positions	-