

# Vessel covers



## Agilent Compatible Vessel Covers

Part No.	Description	OEM Ref
PSCOVERV-00	Low Evaporation Rotary Cover, 0.31" (7.9mm) Cannula Hole Diameter	–
PSCOVERV-01	Clear Cover with Paddle Slot, 0.31" (7.9mm) Cannula Hole Diameter	12-6300
PSCOVERV-02	Amber Cover with Paddle Slot, 0.31" (7.9mm) Cannula Hole Diameter	12-6320
PSCOVERV-03	Clear Solid Cover	12-6330
PSCOVERV-04	Clear Evaporation Cover for Dosage Delivery Module, 0.31" (7.9mm) Cannula Hole Diameter	12-6310
PSCOVERV-HG	Clear Low Evaporation Hinged Cover, 0.31" (7.9mm) Cannula Hole Diameter	–
PSCOVERV-LL	Low Loss Evaporation Cover, Conical Shape, 0.78" (19.8mm) Cannula Hole Diameter	12-0469
PSCOVERV-LLHG	Low Loss Evaporation Hinged Cover, Conical Shape, 0.78" (19.8mm) Cannula Hole Diameter	12-6328
PSCOVERV-LLSB	Stationary Basket Apparatus Low Loss Evaporation Hinged Cover, Conical Shape	–
PSCOVERV-LLVS	TruCenter Low Loss Evaporation Hinged Cover, Conical Shape, 0.78" (19.8mm) Cannula Hole Diameter	12-6327
PSCOVERV-LL200	Mini Vessel Cover, Conical, Low Loss for 200ml vessel	–
PSCOVHGV-SB	Clear Stationary Basket Apparatus Hinged Cover	–
PSCOVERV-6K	Clear Evaporation Cover for VK600/VK6000 Series, 0.31" (7.9mm) Cannula Hole Diameter	12-6200

## Distek Compatible Vessel Covers

Part No.	Description	OEM Ref
PSCOVERD-00	Low Evaporation Rotary Cover, 0.53" (13.5mm) Cannula Hole Diameter	–
PSCOVERD-01	Clear Hinged Cover, 0.53" (13.5mm) Cannula Hole Diameter	3250-0089
PSCOVERD-02	Amber Hinged Cover, 0.53" (13.5mm) Cannula Hole Diameter	–
PSCOVERD-BK	Clear Hinged Cover with Large Centre Hole for Easy Basket Access, 0.53" (13.5mm) Cannula Hole Diameter	3250-0070; 2800-0870
PSCOVERD-ABK	Amber Hinged Cover with Large Centre Hole for Easy Basket Access, 0.53" (13.5mm) Cannula Hole Diameter	–
PSCOVERD-HG	Clear Low Evaporation Hinged Cover, 0.53" (13.5mm) Cannula Hole Diameter	–
PSCOVHDG-SB	Clear Stationary Basket Apparatus Hinged Cover	3250-0081
PSCOVPLG-7D	7th Position Cover Plug	3350-0109
PSCOVPLG-DK	Probe Hole Plug, pk/2	5720-4037
PSCOVCAP-DK	Clear Cap for Basket Access Cover	3250-0019; 2800-0055
PSCOVCAP-ADK	Amber Cap for Basket Access Cover	–
PSHINGE-DK	Plastic Hinge	5560-6011

## Hanson Compatible Vessel Covers

Part No.	Description	OEM Ref
PSCOVERH-00	Low Evaporation Rotary Cover, 0.53" (13.5mm) Cannula Hole Diameter	65-233-170
PSCOVERH-01	Clear Cover with Paddle Slot	72-500-215
PSCOVERH-02	Grey Cover with Paddle Slot	72-800-603
PSCOVERH-03	Amber Cover with Paddle Slot	–
PSCOVERH-04	Clear Solid Cover	–
PSCOVERH-HG	Clear Low Evaporation Hinged Cover, 0.53" (13.5mm) Cannula Hole Diameter	95-590-046
PSCOVERH-LL	Low Loss Evaporation Cover, Conical Shape	–



PSCOVERD-00



PSCOVERD-01



PSCOVERD-02



PSCOVERD-BK



PSCOVERD-ABK



PSCOVERD-HG



PSCOVDHG-SB



PSCOVPLG-7D



PSCOVPLG-DK



PSCOVCAP-DK



PSHINGE-DK



PSCOVERH-00




PSCOVERH-01



PSCOVERH-02



PSCOVERH-03



PSCOVERH-04



PSCOVERH-HG




PSCOVERV-00



PSCOVERV-01



PSCOVERV-02



PSCOVERV-03



PSCOVERV-04



PSCOVERV-HG



PSCOVERV-LL



PSCOVERV-LLHG



PSCOVERV-LLVS



PSCOVERV-6K



PSCOVERV-LL200