

## Serological Pipettes

### Sterilin Plugged Pipettes, Sterile, Individually Wrapped in Paper Peel Packaging

- Printed with product code/lot number and expiry date
- Paper backing allows for easier access
- Paper backing maintains shape for easy disposal of pipette within wrapper following use
- Gamma irradiated to ensure sterility
- Certified non-pyrogenic with an endotoxin level of below 0.25EU/ml
- Batch specific certificate supplied in each carton
- Lot number for complete traceability
- Accuracy of  $\pm 1\%$  for nominal volume
- Clear black printed graduations for maximum clarity
- Descending and ascending graduations to aid pipetting
- Negative graduations for extra capacity
- Colour-coded for ease of identification



### Sterilin Pipettes, Sterile, Individually Wrapped in Paper Peel Packaging

Catalogue Number	Description	Negative Graduations (ml)	Overall Length (mm)	Sterility	Pack Size
RC1640	1ml sterile pipette individually wrapped	-0.2	279	IRR	10 x 100
RC1650	2ml sterile pipette individually wrapped	-0.3	279	IRR	10 x 100
RC1030	5ml sterile pipette individually wrapped	-2.0	275	IRR	4 x 50
RC1040	10ml sterile pipette individually wrapped	-3.0	333	IRR	10 x 50
RC1050	25ml sterile pipette individually wrapped	-10.0	349	IRR	8 x 25
RC1060	50ml sterile pipette individually wrapped	-5.0	344	IRR	2 x 25


### Sterilin Plugged Pipettes, Bulk Packed

- Bulk packed in convenient shelf packs
- Gamma irradiated to ensure sterility
- Certified non-pyrogenic with an endotoxin level of below 0.25EU/ml
- Batch specific certificate supplied in each carton
- Lot number for complete traceability
- Accuracy of  $\pm 1\%$  for nominal volume
- Clear black printed graduations for maximum clarity
- Descending and ascending graduations to aid pipetting
- Negative graduations for extra capacity
- Colour coded for ease of identification



### Sterilin Pipettes, Sterile, Bulk Packed

Catalogue Number	Description	Negative Graduations (ml)	Overall Length (mm)	Sterility	Pack Size
RC1021	1ml sterile pipettes bulk packed	-0.2	279	IRR	20 x 50
RC1022	2ml sterile pipettes bulk packed	-0.3	279	IRR	20 x 50
RC1025	10ml sterile pipettes bulk packed	-3.0	333	IRR	20 x 25
RC1023	5ml sterile pipettes without suction adapter	-1.5	275	IRR	8 x 25



## Sterilin Plugged Pipettes, Individually Wrapped in Plastic Film Packaging

- Gamma irradiated to ensure sterility
- Certified non-pyrogenic with an endotoxin level of below 0.25EU/ml
- Batch specific certificate supplied in each carton
- Lot number for complete traceability
- Accuracy of  $\pm 1\%$  for nominal volume
- Clear black printed graduations for maximum clarity
- Descending and ascending graduations to aid pipetting
- Negative graduations for extra capacity
- Colour-coded for ease of identification



### Sterilin Pipettes, Sterile, Individually Wrapped in Plastic Film Packaging

Catalogue Number	Description	Negative Graduations (ml)	Overall Length (mm)	Sterility	Pack Size
RC1061	1ml sterile pipette individually wrapped	-0.2	279	IRR	10 x 100
RC1062	2ml sterile pipette individually wrapped	-0.3	279	IRR	10 x 100
RC1064	5ml sterile pipette individually wrapped	-2.0	298	IRR	4 x 50
RC1065	10ml sterile pipette individually wrapped	-3.0	333	IRR	10 x 50
RC1066	25ml sterile pipette individually wrapped	-10.0	349	IRR	8 x 25
RC1067	50ml sterile pipette individually wrapped	-5.0	344	IRR	2 x 25
RC1068	5ml sterile pipette without suction adapter	-1.5	275	IRR	4 x 50

