

## PCR Tubes

Designed for optimal heat transfer and sample security, this range of 0.2ml and 0.5ml PCR tubes offers excellent performance and value.

Choose from domed or flat caps. PCR tubes are also available in a range of colours so you can colour co-ordinate samples in your lab to aid workflow.

- Designed to fit all major suppliers' instruments
- 0.2ml or 0.5ml tubes available
- Choice of caps
  - Flat caps for easy sample labelling
  - Domed caps for better heat transfer from thermal cycler lids
- Certified free from RNase, DNase, DNA, DNA Inhibitors and endotoxins



### Individual 0.2ml PCR Tubes

Catalogue Number	Description	Pack Size
LW2330*	0.2ml PCR Tubes, domed caps	1000
LW2170	0.2ml PCR Tubes, flat caps	1000

\*Also available in colours. Please see our website.

### Individual 0.5ml PCR Tubes

Catalogue Number	Description	Pack Size
LW2175*	0.5ml PCR Tubes, flat caps	1000

\*Assorted packs available, add suffix: AS=assorted.


## PCR Tube Strips

Our range of PCR strip tubes has been selected to bring you a comprehensive range for standard and Quantitative PCR.

With all the features of our standard PCR tubes, the strip tube range now includes strips of white tubes with individually attached optically clear caps.

The white Wells maximise reflectance for enhanced signals whilst the optically clear caps enable optimal signal transmission and detection. Optically clear caps are also available on clear strip tubes.

- Some formats can be cut into individual tubes or groups of tubes
- Available without caps for use with cap strips
- Optical, flat or domed caps available
- White tubes for enhanced signals, clear tubes for sample visibility



### 8- and 12-Strip Thin-Walled Tubes and Caps for PCR

Catalogue Number	Description	Pack Size
<b>Strip tubes without caps</b>		
LW2520	0.2ml 8-strip tubes, without caps	125
LW2161	0.2ml 12-strip tubes without caps	80 x 12
<b>Strip tubes with separate strips of caps</b>		
LW2335	0.2ml 8-strip tubes & 8-strip domed caps	125
LW2535	0.2ml 8-strip tubes & 8-strip flat optical caps	125
LW2340	0.2ml 12-strip tubes & 12-strip domed caps	80
<b>Strip tube with attached hinged strip of caps</b>		
LW2500	0.2ml 8-strip tubes & hinged 8-strip domed cap	125
LW2560	0.2ml 8-strip tubes & hinged 8-strip flat optical cap	125
<b>Strip tubes with individually attached caps</b>		
LW2510*	0.2ml 8-strip tubes & individual domed caps	120
LW2505	0.2ml breakable 8-strip tubes & individual flat optical caps	120
LW2570*	0.2ml 8-strip tubes & individual flat optical caps	120
LW2575	0.2ml 8-strip white tubes & individual flat optical caps	120
LW2580	0.1ml 8-strip low profile tubes & individual flat optical caps	10 x 12
<b>Cap Strips - for use with strip tubes above and plates as shown in the compatibility chart</b>		
LW2222	8 Strip Domed Caps, Natural	125
LW2141	0.2ml 8-Strip Caps, Optically Clear	125
LW2190	0.2ml 12-Strip Domed Caps Natural	80
LW2525	0.2ml 12-Strip Flat Optical Caps	80
<b>Cap Strips - for use with plates as shown in the compatibility chart</b>		
LW2905	Domed cap strip, 8, for PCR plates	125
LW2915	Flat cap strip, 8, for PCR plates	125
LW2925	Cap strip, 12, for PCR plates	125

\* Also available in colours. Please see our website.

### Corbett RotorGene qPCR Strip Tubes with Caps

Catalogue Number	Description	Pack Size
LW2590	0.1ml qPCR tubes with caps 4 strip	250
LW2590-BULK	0.1ml qPCR tubes with caps 4 strip	10 x 250