

TEKNOKROMA	PHASE COMPOSITION	AGILENT	SUPELCO	RESTEK	AGILENT	SGE	ALLTECH	QUADREX	PHENOMENEX	MACHEREY-NAGEL	USP CODE
SAPIENS-1MS	100% dimethyl polysiloxane	DB1-MS UI, HP-1 MS UI, VF-1MS	Equity-1	Rxi-1MS	VF-1MS	SOL-GEL-1MS, BPX-1			ZB-1MS	OPTIMA-1MS Accent	
SAPIENS-5MS	95% dimethyl-5% diphenylpolysiloxane	HP-5 MS UI	Equity-5	Rxi-5MS		BP-5			ZB-5MSi	OPTIMA-5MS	
SAPIENS-X5MS	polysiloxane containing p-silphenylene	DB-5MS UI, VF-5MS	SLB-5MS	Rxi-5Sil MS		BPX-5, BP-5MS			ZB-5MS	OPTIMA-5MS Accent	
SAPIENS-WAX.MS	100% polyethylene glycol	VF-WAXMS	SupelcoWAX 10	STABILWAX MS	VF-WAXMS	SolGel-WAX			ZB-WAX	Optima-WAX	G14, G15, G16, G20, G39
SAPIENS-WAX.HT	polyethylene glycol										
TRB-1	100% dimethyl polysiloxane	HP-1, HP101, Ultra-1, DB-1	SPB-1, Equity-1	Rtx-1	CP-Sil 5 CB	BP-1	AT-1	007-1	ZB-1	Optima-1	G1,G2,G38
TRB-1HT		DB-1HT		Stx-1HT			AT-1 HT		ZB-1HT		
TRB-1HT SimDist		DB-1HT SimDist		MXT-1 SimDist	CP-SimDist						
TRB-1MS		HP-1MS, DB-1MS	EQUITY-1	Rxi-1MS, Rtx-1MS	VF-1MS, CP-Sil 5 CB MS		AT-1 MS		ZB-1MS	Optima-1MS	
TRB-SULFUR			SPB-1Sulfur		CP-Select CB for Sulfur						
TRB-2887		DB-2887	Petrocol-2887	Rtx-2887							
TRB.50.2PONA		HP-PONA	Petrocol DH 50.2	Rtx-1 PONA	CP-SIL PONA CB	BP-1 PONA					
TRB-PETROL		DB-PETRO	Petrocol DH	Rtx-1 PONA							
TRB-PETRO.150			Petrocol DH 150								
TRB-5	95% dimethyl-5% diphenyl polysiloxane	HP-5, Ultra-2, DB-5	SPB-5, Equity-5	Rtx-5	CP-SIL 8 CB	BP-5	AT-5	007-2	ZB-5	Optima-5	G27, G36
TRB-5HT		DB-5HT							ZB-5HT		
TRB-5MS		HP-5MS, PAS-5	Equity-5	Rtx-5 MS, Rxi-5MS	VF-5MS, CP-Si 8 CB MS				ZB-5MS	Optima-5MS	
TRB-STEROL			SAC-5								
TRB-5 AMINE			PTA-5	Rtx-5Amine	CP-Si 8 CB MS					Optima-5 Amine	
TRB-5.625	95% dimethyl-5% diphenyl polysiloxane	DB-5.625	PTE-5								
TRB-G27			G27	Rtx-G27							
MTI-5		HP-5MSi	PTE-5	Rxi-5MS							
Meta.X5	95% dimethyl-5% polysilphenylene	HP-5TA, DB-5MS	MDN-5, SLB-5MS	Rxi-5Sil MS	CP-Sil 8 CB Low Bleed/MS VF-5MS	BPX-5	AT-5MS	007-5 MS	ZB-5MS	Optima-5 Accent	
Meta.XLB	Proprietary bonded phase	DB-XLB	MDN 12		VF-XMS				ZB-XLB	OPTIMA XLB	
Meta.WAX 400	100% polyethylene glycol				CP-Carbowax 400						

TEKNOKROMA	PHASE COMPOSITION	AGILENT	SUPELCO	RESTEK	AGILENT	SGE	ALLTECH	QUADREX	PHENOMENEX	MACHEREY-NATHAN	USP CODE
Meta.X5 TRIAZINE	proprietary bonded phase										
Meta.VOC	proprietary bonded phase	DB-502.2, HP-VOC	VOCOL	Rtx-502.2							
TRB-1301,TRB-624,	6% cyanopropylphenyl-94% dimethyl polysiloxane	HP-1301,HP-624	SPB-1301,OVI-43	Rtx-1301,Rtx-624	CP-Select 624 CB	BPX-624	AT-624		ZB-624	Optima-1301	G43
TRB-G43		DB-1301,DB-624								Optima-624	
TRB-14	14% diphenyl-86% dimethyl polysiloxane				CP-Sil 13 CB						
TRB-20	20% diphenyl-80% dimethyl polysiloxane		SPB-20,VOCOL	Rtx-20			AT-20	007-7			G28,G32
TRB-35	35% diphenyl-65% dimethyl polysiloxane	HP-35,DB-35	SPB-35	Rtx-35		BPX-35,BPX-608	AT-35	007-11	ZB-35		G42
TRB-1701	14% cyanopropylphenyl-86% dimethyl polysiloxane	HP-1701,PAS-1701 DB-1701	SPB-1701	Rtx-1701	CP-Sil 19 CB	BP-10	AT-1701	007-1701	ZB-1701	Optima-1701	G46
TRB-225	50% cyanopropylphenyl-50% dimethyl polysiloxane	HP-225, DB-225		Rtx-225	CP-Sil 43 CB	BP-225	AT-225	007-225		Optima-225	G7,G19
TRB-PAG	50% polyethylene-50% poly propylenglycol		PAG								
TRB-FFAP	treated polyethylene glycol for acidic compounds	HP-FFAP,DB-FFAP	Nukol, SP-1000	StabilWAX-DA	CP-WAX 58 CB	BP-21	AT-1000,FFAP	007-FFAP		Permabond FFAP	G25,G35
TRB-50	50% diphenyl-50% dimethyl polysiloxane	HP-50+ DB-17,DB-608	SPB-50, SPB-2250	Rtx-50, Rxi-17	CP-Sil 24 CB		AT-50	007-17	ZB-50	Optima-17	G3
TRB-50HT	50% diphenyl-50% dimethyl polysiloxane	DB17HT		Rtx-65	TAB-CB			007-65HT			G17
TRB-F50	50%trifluoropropylmethyl polysiloxane	DB-210,DB-200		Rtx-200			AT-210	007-210		Optima-210	G6
TRB-WAX	100% polyethylene glycol	HP-20M,HP-imoWAX DB-WAX,DB-WAXetr	SupelcoWAX-10, CarboWAX 20M	StabilWAX	CP-WAX 52 CB	BP-20	AT-WAX	007-CW	ZB-WAX	Permabond CW 20M	G14,G15,G16
SupraWAX-280	100% polyethylene glycol		SupelcoWAX-10								
TR-WAX.DB	treated polyethylene glycol for basic compounds	CAM,HP-BasicWAX	CarboWAX-Amine	StabilWAX-DB	CP-WAX 51 CB						
Meta.WAX	100% polyethylene glycol	HP-WAX,DB-WAX			CP-WAX 57 CB						
TRB-WAX.Omega	100% polyethylene glycol		OmegaWAX	FameWAX							
TR-CN100	100% biscyanopropyl polysiloxane		SP-2340	Rt-2340	CP-Sil 88	BPX-70, BPX-90					G48
TR-CRESOL	not bonded phase				CP-Cresol						
TRB-17	50% diphenyl-50% dimethyl polysiloxane	HP-17									G3
Meta.VOC	bonded phase	DB-502.2,HP-VOC	Vocol	Rtx-502.2							
Meta.BLOOD 1	bonded phase	DB -ALC1		Rtx-BAC1							
Meta.BLOOD 2	proprietary bonded phase	DB-ALC2		Rtx-BAC2							
TRB-608	proprietary bonded phase	HP-608	SPB-608			BP-608					
TR-TCEP	1,2,3-tris(cyanoethoxy) propane		TCEP	Rt-TCEP	CP-TCEP						
MetAmine-VOL				Rtx-Volatile Amine	CP-Volamine						
TRB-BIODIESEL	proprietary bonded phase	DB-5HT		Rtx-Biodiesel	Select Biodiesel			400-5HT	ZB-5HT		