

Capsule Sinkers



Sinker Type	Part Number	Style	Material	Coating	Inside Diameter		Inside Length		Wire Diameter		Max Capsule Size
					(inches)	(mm)	(inches)	(mm)	(inches)	(mm)	
BASKET	PSCUSBK-JP	Basket, JP	316 SS	-	0.470	11.9	1.000	25.4	0.040	1.00	00
	PSBSKoo8-01	Basket, 8 Mesh	316 SS	-	0.620	15.7	1.060	26.9	0.028	0.71	000
	PSBSKoo8-JP	Basket, 8 Mesh	316 SS	-	0.510	12.9	0.900	22.9	0.028	0.71	0
	PSSNKoo8-VK	Basket, 8 Mesh	316 SS	-	0.795	20.2	0.500	12.7	0.028	0.71	1
	PSSNKo10-VK	Basket, 10 Mesh	316 SS	-	0.795	20.2	0.500	12.7	0.025	0.64	1
	PSSNKo20-VK	Basket, 20 Mesh	316 SS	-	0.795	20.2	0.500	12.7	0.016	0.41	1
	PSSNKo40-VK	Basket, 40 Mesh	316 SS	-	0.788	20.2	0.500	12.7	0.010	0.25	1
	PSSNKo80-VK	Basket, 80 Mesh	316 SS	-	0.795	20.2	0.500	12.7	0.007	0.18	1
PRONGED	PSCAPWHT-VK	3 Prong	Nylon / PP	-	0.375	9.5	1.000	25.4	-	-	1 (min=3)
	PSMAGWHT-VK	3 Prong	Nylon / PP	-	0.375	9.5	1.000	25.4	-	-	1 (min=3)
	PSCAPWHT-SV	3 Prong	316 SS	-	0.375	7.0	1.000	25.4	-	-	1 (min=3)
O-RING	PSCAPWHT-599	O-Ring	316 SS	-	0.230	5.8	0.900	22.9	-	-	3
SPIRAL	PSCAPWHT-TR	Butterfly	316 SS	-	0.211	5.4	0.375	9.5	0.030	0.76	-
	PSCAPWHT-PC	Strip Film	316 SS	-	0.425	10.8	0.425	10.8	0.037	0.94	-
	PSCAPWST-15	Helix	316 SS	-	0.197	5.0	0.610	15.5	0.016	0.41	5
	PSCAPWHT-04	Helix	Music Wire	PTFE	0.385	9.8	0.840	21.3	0.043	1.09	0
	PSCAPWHT-45	Helix	316 SS	-	0.385	9.8	0.840	21.3	0.047	1.19	0
	PSCAPWHT-LG	Helix	Music Wire	PTFE	0.410	10.4	1.100	27.9	0.043	1.09	000
	PSCAPWHT-XL	Helix	Music Wire	PTFE	0.465	11.8	1.150	29.2	0.043	1.09	000
	PSCAPWHT-XLS	Helix	316 SS	-	0.465	11.8	1.150	29.2	0.043	1.09	000
	PSCAPWHT-XLSS	Helix	316 SS	-	0.465	18.8	0.885	22.5	0.043	1.09	000
	PSCAPWST-18	Helix	316 SS	-	0.236	6.0	0.709	18	0.016	0.41	Special
	PSCAPWST-19	Helix	316 SS	-	0.276	7.0	0.760	19.3	0.016	0.41	3
	PSCAPWHT-25	Helix	316 SS	-	0.370	9.4	0.840	21.3	0.043	1.09	2
	PSCAPWHT-02	Helix	Music Wire	PTFE	0.370	9.4	0.840	21.3	0.047	1.19	0
	PSCAPWHT-XS	Helix	316 SS	-	0.575	14.6	0.885	22.5	0.047	1.19	0
	PSCAPWHT-23	Helix	316 SS	-	0.315	8.0	0.905	23	0.016	0.41	Special
PSCAPWHT-31	Helix	316 SS	-	0.433	11.0	1.220	31	0.031	0.78	0	
WIRE	PSWIRESK-VK	Sinker Wire 50'	316 SS	-	-	-	-	0.032	0.81	000	

CAPSULES				
Capsule Size	Length		Diameter	
	(inches)	(mm)	(inches)	(mm)
000	1.029	26.1	0.390	9.91
00E	1.012	25.7	0.336	8.53
00	0.917	23.3	0.336	8.53
0E	0.909	23.1	0.301	7.65
0	0.854	21.7	0.301	7.65
1	0.765	19.4	0.272	6.91
2	0.709	18.0	0.250	6.35
3	0.626	15.9	0.229	5.83
4	0.563	14.3	0.209	5.32
5	0.437	11.1	0.193	4.91

